





DISS Conversion Outlets - Puritan Series A

M3524A Oxygen | M3522 Vacuum | M3516 Medical Air | M3504 Nitrous Oxide

Puritan Series 73 Outlet



Class 1 Conversion Outlet




Tools



Phillips Screwdriver



5/8" Hex Wrench



Flat Blade Screwdriver



5/16" Hex Wrench



5/32" Hex Key Wrench

Refer to reverse side for detailed instructions.



DISS Conversion Outlets - Puritan Series A

M3524A Oxygen | M3522 Vacuum | M3516 Medical Air | M3504 Nitrous Oxide

Puritan Series A NPT Outlet Disassembly

<p>1</p> <p>CLOSE THE AREA ZONE VALVE.</p> <p>Warning: Do not close zone valve until all personnel have been advised of the intended service and all systems requiring medical gas are being supplied from portable or alternate supplies.</p>	<p>2</p> <p>Remove the mounting screws, coverplate, and trim plate.</p>	<p>3</p> <p>HOLDING SCREW</p> <p>Loosen the holding screw on the valve block using a 5/16" hex wrench.</p>	<p>4</p> <p>Turn the quick-connect counterclockwise to unscrew from the valve block.</p>
<p>5</p> <p>Remove the quick-connect and the spring and poppet.</p>	<p>6</p> <p>Remove the screws and indexing collar from the valve block.</p>	<p>7</p> <p>Remove the screws and flange from the valve block.</p>	<p>8</p> <p>Using a 5/8" hex wrench, turn the connector counterclockwise to remove the connector and poppet.</p>

Class 1 Conversion Outlet Installation Instructions

<p>1</p> <p>Clean all loose chips and debris from the valve block.</p>	<p>2</p> <p>SEAL RINGS</p> <p>Insert the conversion outlet into the valve block cavity being careful not to damage the seal rings on the outlet.</p>	<p>3</p> <p>APPLY HERE</p> <p>Apply Loctite® adhesive to the threads on the cap screws and use a 5/32" hex key wrench to attach the outlet to the valve block.</p>	<p>4</p> <p>Test outlet in accordance with CSA Pamphlet Z305.1.</p>
<p>5</p> <p>Align sub-mounting plate so that the gas name can be read from the front. Attach with mounting screws.</p>	<p>6</p> <p>INDEXING PIN</p> <p>INDEXING HOLE</p> <p>Align the coverplate so that the indexing pin slides through the indexing hole on the sub-mounting plate. Attach with mounting screws.</p>	<p>7</p> <p>Conversion complete. Attach mating equipment.</p>	