

Hose Order Guide



Gas:

O2 = Oxygen

AIR = Medical Air

VAC = Medical Vacuum

N2O = Nitrous Oxide

N2 = Nitrogen

CO2 = Carbon Dioxide

AGSS = Anaesthetic Gas Scavenging System

WAGD = Waste Anesthetic Gas Disposal
EVAC = Laser Smoke Evacuation


IA = Instrument Air

Inlet Connections:

= Female Wrench Tight DISS FΗ = Female Hand-Tight DISS = Male DISS Check Valve M MB = Male DISS Body Adaptor PF = Puritan Female PM = Puritan Male OF = Ohio Female OM = Ohio Male XF = Oxequip Female XM = Oxequip Male MF = Medstar Female MM = Medstar Male CF = Chemtron Female = Chemtron Male CM SF = Schrader Female SM = Schrader Male

Outlet Connections:

F = Female Wrench Tight DISS FΗ = Female Hand-Tight DISS M = Male DISS Check Valve MB = Male DISS Body Adaptor = Puritan Female = Puritan Male PM= Ohio Female OF OM = Ohio Male XF = Oxequip Female XM= Oxequip Male = Medstar Female MF MM = Medstar Male CF = Chemtron Female CM = Chemtron Male SF = Schrader Female SM = Schrader Male



Example-AIR-FHM-3

